

# SynDaver Labs

## SynDaver Surgical Canine

The SynDaver Surgical Canine is a revolutionary high fidelity surgical trainer. The model provides an unparalleled platform for repeatable surgical procedures. Featuring customizable pathology integrated with complete vasculature, the system affords surgeons an incomparable experience in surgical situations commonly faced by veterinarians.

Capabilities of the SynDaver Surgical Canine include:

### Abdominal Procedures:

- Gastrotomy with foreign body removal
- Gastropexy
- Splenectomy
- Liver lobe biopsy
- Partial liver lobectomy
- Enterotomy with foreign body removal
- Intestinal resection with anastomosis
- Small intestinal biopsy
- Large intestinal biopsy
- OVH
- Cystotomy with stone removal


### Other Surgical Procedures:

- Tracheostomy
- Mammary mass resection
- Cutaneous suturing
- Muscular suturing
- Surgical hemostasis

### Technical Procedures:

- Cephalic venipuncture
- Cephalic vein catheterization
- External jugular venipuncture
- External jugular catheterization
- Orotracheal intubation
- Nasoesophageal tube placement
- Stomach tube placement
- Thoracocentesis
- Cystocentesis
- Chest tube placement
- Urinary catheterization
- Upper GI endoscopy
- Lower GI endoscopy
- Cystoscopy

Scan the code  
for videos of the  
SynDaver Surgical  
Canine in action:



# SynDaver Labs

8506 Benjamin Road  
Tampa, FL 33634  
813.600.5530  
syndaver.com

