

The Adult Cric Trainer is the world's most realistic surgical training platform for cricothyrotomy, nasogastric intubation, and retrograde intubation. This model allows students to practice and repeat technique on a high quality, live-tissue replacement platform in a biohazard-free environment.

Repetitive use will strengthen the ability and confidence of all team members who perform or assist in implementing surgical airways. Anatomical features not shown include oral cavity and nasal passages communicating with the lower airway.

Relevant Skills

Surgical and needle cricothyrotomy, nasal intubation, retrograde airway, palpation, cannulation, administration and removal of sutures and staples, surgical cutdown and application of adhesives and bandages.

Equipment Compatibility

Surgical airway devices, autosuturing and autostapling devices, laser scalpels, bipolar and monopolar devices, electrocautery devices, harmonic blades and all imaging equipment.


Tissue Features

SynTissue synthetic human tissues made from salt, water and fiber, which feature the world's most realistic tactility. SynTissue synthetic human tissues match the acoustical characteristic of real human tissue.

Cric Replacement Tissues

Each Adult Cric Trainer includes 20 replacement cric skin tissues and membranes.

Additional replacements are available in packs of 30 skin replacement tissues and 30 cricoid membranes. [Part number: 160710]



Included Components

- Base Board
- Muscular Phantom
- Skin Overlay
- Replaceable Cric Carriages (x2)
- Replacement Tissue Sets (x20)
- Storage Case
- Care Kit