
# SynDaver Labs

Phone: 813.600.5530 Fax: 813.600.3235 info@syndaver.com www.syndaver.com

SynDaver Liver

130610



Our SynDaver liver models are by far the most realistic synthetic organs of their type available anywhere in the world. The structural design is based on an amalgam of CT and MRI images from actual patients and the synthetic tissues employed in construction have been validated against the mechanical, physicochemical, thermal and dielectric properties of living tissue.

This organ is one of many available separately from our SynDaver Synthetic Human product line. If you are looking for a custom version of this organ or one which is not listed please call us. We can build organs from CT and MRI images, CAD designs or simple drawings.

# **Organ Features**

Skinned outer shell with lobed structure. This organ is also available with the primary arterial and venous trunks, complete biliary system and a variety of pathologies.

# **Options**

Available with biliary system options. If you require additional customization please call and ask to speak to one of our technical sales representatives.

# **Typical Uses**

These organs are used in the SynDaver Synthetic Human product line. They are also incorporated into complex model systems for liver transplant training and biliary stent testing.

### **Imaging Equipment**

Compatible with all known imaging equipment including MRI, CT, fluoroscopy and ultrasound.

### **Surgical Equipment**

Compatible with all known surgical devices including lasers, RF ablation, bipolar, monopolar and harmonic devices.

### **Extraordinary Features**

SynTissue synthetic human tissues made from salt, water and fiber—which feature the world's most realistic tactility. SynTissue synthetic human tissues match the acoustical characteristic of real human tissue.

All of our products are made in the USA.