

SynDaver's® Airway Trainers are realistic medical training platforms ideal for teaching the techniques associated with tracheal intubation.

With these trainers, students will be able to learn and master surgical techniques on biohazard-free material that looks, feels and behaves like live human tissue.

These models include realistic oral cavity with a hard and soft palate, tongue, uvula, epiglottis and vocal cords. The soft neck with cricocartilage enables users to perform Sellick's maneuver for a better view of the larynx and/or to reduce gastric reflux.

Relevant Skills

- Intubation
- Airway Management Exercises

Components Included

- Torso with age-accurate anatomy
- Storage Case
- Care Kit

Tissue Features

SynTissue® is made from water, fibers and salts, boasting the world's most lifelike and realistic tactility next to living biological tissues. SynTissue synthetic tissues match the physical and mechanical properties of live tissue.



For More Information About SynDaver Products
813.600.5530 • SynDaver.com