



Relevant Skills

Users can practice surgical and needle cricothyrotomy, nasal intubation, retrograde airway, palpation, cannulation, application and removal of sutures and staples, surgical cutdown and application of adhesives and bandages.

Components Included

Plastic base, muscular form, two membrane carriages with hyoid and cricoid, skin overlay, 20 replacement tissues (skin and cricoid membrane) and a durable Pelican Case.

Options Available

Choice of skin tone and replacement tissue sets (skin and cric membrane) is available.

Equipment Compatibility

SynDaver products are compatible with surgical airway devices, autosuturing and autostapling devices, laser scalpels, electrocautery devices, bipolar and monopolar devices, harmonic blades and all known imaging equipment.

Our SynAtomy® Adult Surgical Cric Trainer is the world's most realistic surgical training platform for cricothyrotomy. In addition, it is compatible with nasogastric intubation and retrograde intubation. This model enables students to practice and repeat and improve their techniques on a high-quality, live-tissue replacement platform in a biohazard-free environment.

Repetitive use strengthens skills and confidence of all team members who perform or assist in implementing surgical airways. Typical users include emergency medical technicians, flight nurses, combat medics, ICU nurses and nurse practitioners. Anatomical features include oral cavity and nasal passages communicating with the lower airway, chin, clavicle, hyoid bone, thyroid cartilage, cricoid cartilage and cricoid membrane.

For More Information About SynDaver Products
813.600.5530 • SynDaver.com

SynDaver
Beyond Human