





www.SynDaver.com



Canine Airway

SynDaver's ultra-high-fidelity canine intubation and upper airway trainer. This model contains highly realistic oral and upper airway anatomy and is by far the most realistic and elaborate canine model that has ever been constructed. This model will allow for an unparalleled intubation experience.



Relevant Skills

The airway offers the ability for respiration if hooked up to a respirator. Cuff inflation is accessible as well as pressure checks and leak detection. Inflatable lungs are available on the deluxe model. Both models support endotracheal intubation



Included Components

Larynx, epiglottis, glottis, and tracheal opening. The deluxe model contains a vascularized neck skin to allow for IV.



Available Options

If you require custom tissues, dimensions, or modified mechanical properties, please call and ask to speak to one of our technical sales representatives.



Packaging

Product is shipped vacuum packed with storage case.

