

Adult Cricothyrotomy Trainer

SynDaver Adult Surgical Cric Trainer is the most realistic platform available for cricothyrotomy training. This model allows students to practice, repeat, and refine their techniques on a high-fidelity, live-tissue replacement platform in a biohazard-free environment.

Repetitive practice on this model builds the skills and confidence needed by all team members involved in surgical airway procedures. Ideal for emergency medical technicians, flight nurses, combat medics, ICU nurses, and nurse practitioners, this trainer is designed to simulate key anatomical structures, including the lower airway, chin, clavicle, hyoid bone, thyroid cartilage, cricoid cartilage, and cricoid membrane.



Relevant Skills



Users can practice surgical and needle cricothyrotomy, retrograde airway, palpation, cannulation, application and removal of sutures and staples, surgical cutdown and application of adhesives and bandages.

Included Components



Plastic base, muscular form, two membrane carriages with hyoid and cricoid, bifurcated trachea with left and right simulated lung bags, skin overlay, simulated blood, and 30 replacement tissues (skin and cricoid membrane).



Available Options

Replacement tissue sets (skin and cric membrane) and blood analog are available.



Packaging

Product is shipped with accessories in the storage container.