



SynDaver Storage Tank

When not in use, the SynDaver model is kept fully submerged in the table's built-in water chamber to maintain the integrity and longevity of the SynTissue®. When students are ready to interact with the model, one person removes the backboard while two others lift the unit out of the water. The backboard is then slid back into place, allowing the retained water to drain into the storage chamber below.



Construction Specifications

- 2 Piece Powder Coated Aluminum Frame
- High Molecular Weight Polyethylene (HMPE) Tub
- Heliarc welded seams and joints ground and polished to a smooth finish
- Stainless steel trim and hinged hood
- Integrated High Density Polyethylene (HDPE) Backboard
- 6-inch casters, all with brake mechanisms

Key Measurements

Overall Dimensions (L x W x H): 76" x 32" x 36"

Dry Weight: 180 lbs.

Wet Weight (With SynDaver Model Stored): 595 lbs.

Water Capacity: 38 Gallons

Shipping Configuration and Weight:

Box 1 (Base Frame: 44"L x 32.75"W x 16.25"H; 44 lbs.

Box 2 (Upper Frame / Tub): 74.75"L x 34.75"W x 17.75"H; 172 lbs.

