

SynAtomy Deluxe Suturing Kit

160320



Our SynAtomy Deluxe Suturing Kit is a great training platform for both advanced students and experienced professionals seeking to hone the skills associated with more advanced suturing and surgical procedures.

This complete kit includes a great instruction manual, instrument-quality suturing tools, industry-best SynTissue synthetic human tissues, surgical practice board, and a durable travel case.

Relevant Skills: Placement and removal of sutures and staples, anastomosis, laceration repair, cannulation and application of adhesives.

Equipment Compatibility: Laser scalpels, electrocautery and RF ablation devices, harmonic blades, monopolar and bipolar devices, plasma knives, ultrasound equipment and all known imaging equipment.

Included Tools: Surgical practice board, storage case, instruction manual (Suture and Surgical Hemostasis, by Rebecca Pieknick), surgical stapler, staple removal tool, instrument-quality suturing tools (needle driver, scissors, and tweezers), scalpel with 10 replacement blades and 20 assorted sutures.



Included Tissues:

Large basic suture pad, knot tying pad, double layered bowel, two anastomosis vessels, small abdominal suture pad, small muscular suture pad, abdominal aorta and a coronary artery.

Available Options:

If you require custom pathologies, tissues, tools or other changes please call and ask to speak to one of our technical sales representatives.

Extraordinary Features:

SynTissue synthetic human tissues made from salt, water and fiber - which feature the world's most realistic tactility. SynTissue synthetic human tissues match the acoustic characteristic of real human tissue.

All of our products are made in the USA.